

Kat.-Nr./Cat.-No.: C-000.32.549

C-000.32.406

C-000.32.516

C-000.32.405

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.039	VERBINDUNGSELEMENT FASTENERS
02	Z-300.00.040	HALTERUNG KPL. 1 HOLDER 1
03	Z-300.00.025	GEWINDEHÜLSE THREADED STUD
04	Z-300.00.238	DRUCKSCHEIBE THRUST WASHER
05	Z-300.00.346	DREHGRIFF ADJUSTING KNOB
06	Z-300.00.022	HALTERUNG KPL. 2 HOLDER 2
08	C-000.32.403	S-GUARD FÜR KOPFBAND PROF. L S-GUARD FOR HEADBAND PROF.L.
09	C-000.32.401	S-GUARD FÜR KOPFB. LIGHTWEIGHT S-GUARD FOR HEADBAND LW

