

941

**BETA NT Ladegriff AV**

**BETA NT rechargeable handles AV**

Kat.-Nr./Cat.-No.: X-007.99.395

X-007.99.411

X-002.99.411

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.088	GRIFFKOPF AV HANDLE HEAD AV
02	Z-300.00.276	GRIFFHÜLSE OUTER HANDLE
03	X-002.99.382	NIMH-LADEBATTERIE 3,5 V NIMH RECHARG.BATTERY 3,5 V
04	X-007.99.381	LI ION-LADEBATTERIE 3,5V RECH.BATTERY LI ION 3,5V
05	X-002.99.412	BODENEINH.F.BETA-NT-LADEGR. BOTTOM INSERT,BETA NT HANDLE
06	X-007.99.383	LADEBATTERIE LI-ION L RECH. BATTERY LI-ION L
07	X-002.99.396	BODENEINHEIT BETA L BOTTOM INSERT BETA L

