

Kat.-Nr./Cat.-No.:	F-001.22.800 F-008.22.806	F-002.22.806	F-001.22.806	F-008.22.800
--------------------	------------------------------	--------------	--------------	--------------

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	F-000.22.804	GRIFFHUELSE/KL.LARYN.GR. OUTER HANDLE/COMP.LARYN.HA.
02	F-001.22.801	EINSCHUB F.FO/NT-LARYNG.-GR. INSERT SMALL F.O. BATT.HANDLE
03	F-002.22.801	EINSCH. KL.F.O. LARYN. LADEGR. INSERT SMALL F.O.RECH.HANDLE
04	F-008.22.801	LED BATT.EINSCH.F.KL.F.O.LA.GR LED BATT.INS.F.SMALL F.O.LA.HA
05	X-001.88.035	XHL XENON HALOGEN LAMPE 2,5V XHL XENON HALOGEN BULB 2,5V
06	X-002.88.044	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V
07	X-001.99.306	LADEBATTERIE LED M2Z 2,5V RECHARGEABLE BATT. LED M2Z
08	X-002.99.106	NIMH LADEBATTERIE M 3 Z NIMH BATTERY M 3 Z
09	X-001.99.305	NIMH LADEBATTERIE M2Z 2,5V NIMH RECHARGEABLE BATT. M2Z
10	F-000.22.852	BODENEINH. F.KL.F.O.LA.LADE.GR BOTT.CAP F.SMALL FO LED LAR.HA
11	F-000.22.851	BODENKAPPE F.KL.F.O.LA.BATT.GR BOTT.CAP F.SMALL FO LED LAR.HA

