

949

BETA F.O. Laryngoskopgriff

BETA F.O. Laryngoskop handles

Kat.-Nr./Cat.-No.: F-001.22.117 F-002.22.412 F-007.22.412

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.357	GRIFF MODUL 2 HANDLE HEAD MODULE 2
04	Z-300.00.337	SCHALTBUCHSE SWITCH JACK
05	X-001.88.035	XHL XENON HALOGEN LAMPE 2,5V XHL XENON HALOGEN BULB 2,5V
06	X-002.88.044	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V
07	Z-300.00.103	LAMPENEINSATZ 2 BULBHOLDER 2
08	X-007.99.381	LI ION-LADEBATTERIE 3,5V RECH.BATTERY LI ION 3,5V
09	X-002.99.382	NIMH-LADEBATTERIE 3,5 V NIMH RECHARG.BATTERY 3,5 V
10	X-002.99.412	BODENEINH.F.BETA-NT-LADEGR. BOTTOM INSERT,BETA NT HANDLE
11	X-000.99.119	BODENEINHEIT F.BETA BATT.GRIFF BOTTOM CAP F.BETA BATT. HANDLE

