

Kat.-Nr./Cat.-No.: B-001.11.500 B-002.11.500 B-008.11.500

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.199	LINSE LENS
02	Z-300.00.230	DICHTRING SEALING RING
03	Z-300.00.143	BEFESTIGUNGSSCHRAUBE MIT O-RING FIXING SCREW WITH O-RING
05	Z-300.00.229	LINSENDECKEL LENS MOUNT
09	X-001.88.077	XHL XENON HALOGEN LAMPE 2,5V XHL XENON HALOGEN BULB 2,5V
10	X-002.88.078	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V
11	Z-300.00.387	LINSENDECKEL KPL. M. BEFESTIGUNGSSCHRAUBE LENS COVER COMPLETE WITH FIXING SCREW

