

**HEINE BETA® 200 LED Ophthalmoskop**

**HEINE BETA® 200 LED Ophthalmoscope**

Kat.-Nr./Cat.-No.: C-001.30.100 C-002.30.100 C-001.30.102 C-002.30.102

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
18	X-001.88.069	XHL XENON HALOGEN LAMPE 2,5V XHL XENON HALOGEN BULB 2,5V
19	X-002.88.070	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V

