

Kat.-Nr./Cat.-No.: C-000.33.020

C-000.33.327

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.366	STIRNPOLSTER FRONT PAD
08	Z-300.00.025	GEWINDEHÜLSE THREADED STUD
09	Z-300.00.238	DRUCKSCHEIBE THRUST WASHER
10	Z-300.00.239	DREHGRIF FIXATION KNOB
18	Z-300.00.367	KOPFPOLSTER OVERHEAD PAD
19	Z-300.00.215	VERSCHLUßKAPPE COVER
26	Z-300.00.262	KABELHALTER CORD RETAINER
28	Z-300.00.242	KOPFPOLSTER FOAM PAD

