

Kat.-Nr./Cat.-No.: C-001.35.010 C-002.35.010 C-001.35.015 C-002.35.015

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.265	SCHEIBE COVER
02	Z-300.00.258	ZYLINDERSCHRAUBE M2X6 CYLINDER SCREW M2X6
27	Z-300.00.264	STIRNSTÜTZE BROW REST
28	Z-300.00.266	STIRNSTÜTZENPOLSTER BROW REST CUSHION
29	X-001.88.077	XHL XENON HALOGEN LAMPE 2,5V XHL XENON HALOGEN BULB 2,5V
30	X-002.88.078	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V

