

Kat.-Nr./Cat.-No.: **C-004.33.500**

C-004.33.503

C-008.33.502

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.396	INSTRUMENTENKABEL KPL. INSTRUMENT CABLE ASSEMBLY
09	Z-300.00.009	SCHUTZDECKEL PROTECTIVE COVER
16	Z-300.00.214	OKULARMUSCHEL EYEPIECE RUBBER
19	Z-300.00.237	KORREKTURLINSE +2D LENS 2 DIOPTERS
20	Z-300.00.256	LINSE 0 DIOPTRIEN LENS 0 DIOPTRIEN
21	Z-300.00.236	OKULARTRICHTER EYEPIECE FUNNEL
22	X-004.88.111	XHL XENON HALOGEN LAMPE 6V XHL XENON-HALOGEN-HALOGENLAMPE 6V
23	X-008.87.200	LED MODUL FUER OMEGA 500 LED MODULE FOR OMEGA 500

