

Kat.-Nr./Cat.-No.: B-001.11.492

B-002.11.492

Nr. Pos	Bestell-Nr. Order-No.	Bezeichnung Description
01	Z-300.00.197	TRÄGER KPL. SPECULUM HOLDER
02	Z-300.00.048	KUGEL BALL
03	Z-300.00.001	DRUCKFEDER SPRING
04	Z-300.00.002	KORDELMUTTER KNURLED CONTROL SLEEVE
05	X-001.88.037	XHL XENON HALOGEN LAMPE 2,5V XHL XENON HALOGEN BULB 2,5V
05	X-002.88.049	XHL XENON HALOGEN LAMPE 3,5V XHL XENON HALOGEN BULB 3,5V
07	Z-300.00.198	LINSENTRÄGER KPL. LENS HOLDER
08	Z-300.00.216	DRUCKFEDER SPRING
09	Z-300.00.361	ANSATZSCHRAUBE SCREW
13	B-000.11.212	BLLENDE F.DIAGN.OTOSKOP LAMPE SLEEVE F.DIAGN.OTOSC. BULB

